



**Acuity Size**

0 20/400  
 9 20/200  
 1 20/100  
 8 20/80  
 7 20/70  
 6 20/60  
 5 20/50  
 4 20/40  
 3 20/30  
 2 20/20  
 . Reset

**Controls**

H Hide/unhide  
 L Single line  
 U Multi line  
 V Video  
 Space Randomize  
 R, \* Temp. Random  
  
 F4 Random on/off  
 F1 Help Key  
 O One Letter

**Acuity type**

S Snellen  
 I ETDRs  
 P Pictures  
 E Tumbling E  
 C Tumbling C  
 N Numbers  
 T HOTV

**Special test**

W Worth dot  
 F2 Fixations disparity  
 A Astigmatic Display  
 J Reduced Contrast  
 G Red/Green  
 K Crowding Bar

Page Down, - Decrease acuity 1 line  
 Page up - Increase acuity 1 line

**Acuity Size**

<b>0</b>	Acuity 20/400
<b>9</b>	Acuity 20/200
<b>1</b>	Acuity 20/100
<b>8</b>	Acuity 20/80
<b>7</b>	Acuity 20/70
<b>6</b>	Acuity 20/60
<b>5</b>	Acuity 20/50
<b>4</b>	Acuity 20/40
<b>3</b>	Acuity 20/30
<b>2</b>	Acuity 20/20
<b>PgDn, -, Up</b>	Decrease acuity 1 line
<b>PgUp, +, Down</b>	Increase acuity 1 line
<b>Dot</b>	Reset acuity

**Acuity type / special test**

<b>S</b>	Snellen
<b>I</b>	ETDRs
<b>P</b>	Picture
<b>E</b>	Tumbling E
<b>C</b>	Tumbling C
<b>N</b>	Number
<b>T</b>	HOTV
<b>W</b>	Worth dot
<b>F2</b>	Fixation disparity
<b>A</b>	Astigmatic display
<b>J</b>	Reduced contrast
<b>G</b>	Red/Green overlay
<b>K</b>	Crowding bars

**Controls**

<b>H</b>	Hide/unhide
<b>O</b>	One letter
<b>L</b>	Single line
<b>U</b>	Multi line
<b>Space</b>	Randomize
<b>R, *</b>	Temp. random
<b>F4</b>	Random ON/OFF
<b>F1</b>	Help key
<b>V</b>	Video
<b>Scroll</b>	Force Video